

## Motivation

- Automated tracking of Parkinson's disease (PD) motor symptoms during daily activities requires robust and clinically-relevant sensor metrics<sup>1</sup>.
- Bradykinesia, a major motor symptom of PD, is activity-dependent and requires advanced tools for assessment that include gait.
- Bradykinesia can manifest during gait as a reduction in arm swing, leg velocity, range of movement, and heel-toe dynamics, which are difficult to capture with conventional sensors such as accelerometers.
- Recent advances in wearable inertial (IMU) sensors may provide an enhanced means of recording joint kinematic of bradykinesia gait disorders.

## Objective

To evaluate the ability of IMU sensor-based metrics to identify gait impairments associated with body bradykinesia in patients with Parkinson's disease during unscripted activities of daily living.

## Approach

- Body bradykinesia assessment during gait was based on impairments listed in Item 31 of UPDRS.
- An automated classifier (neural network) was designed to isolate gait activity from other activities based on leg sensor data.
- Five gait metrics were analyzed using angular velocity (gyroscope) and compared during presence and absence of bradykinesia.

## Data Collection


### Subject Population

PD Subjects Demographics	
Number	n = 6
Age (y)	57.5 ± 12.5
Male/Female	4/2
Disease Duration(y)	8.6 ± 5.4
Total Data	1000 min
Bradykinesia Prevalence* (%)	58.7
Hoehn-Yahr (ON)	II-III,

\* % of Total Data w/ bradykinesia

### Data Acquisition Protocol

- Data were acquired from 2 wireless IMU sensors (Trigno IM - Delsys Inc) from forearm and shank – see Figure
- Data were recorded continuously for 3 hours in a simulated home setting during unscripted activities.
- Video recordings were acquired and annotated by movement disorder experts to identify activity type and presence/absence of body bradykinesia (based on Item 31 of UPDRS)






## Methods

- Gait impairment metrics were calculated on the basis of angular velocity magnitude and range of movement from gyroscope data. [Refer to Analysis section]
- Statistical comparisons (Mann-Whitney-Wilcoxon test) were computed to test the discriminability of gait impairment metrics for Bradykinesia and Non-Bradykinesia portions of the gait data. [Refer to Results section]


## Analysis

- The numbers [1-5] in the figure refer to the five sensor-based gait metrics that were derived to quantify gait impairments associated with bradykinesia.
- Raw gyroscope data from upper and lower limb are compared with raw Accelerometer data for the same IMU recording.
- These comparisons illustrate the greater precision in identifying data features from the angular velocity plots compared to accelerometer plots.



## Results

- The numbers [1-5] identify the bradykinesia gait impairments described in the Analysis. Corresponding raw angular velocity for instances with and without bradykinesia are analyzed.
- The results demonstrate that IMU sensors provide objective gait metrics with significant discriminability for bradykinesia.



## Conclusion

- The work demonstrates the viability to develop robust and clinically-relevant metrics for improved detection of gait abnormalities in PD.
- A subsequent study [see Adjacent Poster: Roy et al. "Autonomous Tracking of Body Bradykinesia..."] utilized these impairment metrics as features to train a neural network classifier to detect body bradykinesia and achieved <5% error during unscripted activities in a population of n=16 PD subjects.

## Acknowledgments

Research supported in part by the De Luca Foundation and by the National Institute of Neurological Disorders and Stroke of the National Institutes of Health under award R44NS083098.

References <sup>1</sup>Movement Disorders, Vol 31, No 9 2016.